

UVR-M and UVR-Mi, UV Air Recirculators Test Report



UVR-M and UVR-Mi, UV air recirculators Test Report

Author: Dr. biol. Vasily Bankovsky

UV air recirculators **UVR-M** and **UVR-Mi**, produced by Biosan, are equipped with bactericidal UV lamps (Philips) and are used for air disinfection in research laboratories, hospitals and veterinary clinics.

To show the efficiency of UV air recirculators **UVR-M** and **UVR-Mi**, we examined UV intensity in Philips 25W bactericidal UV lamps and an impact of UV radiation on various types of microorganisms.

General information


Theory: UV radiation affects the viability of microorganisms by causing photochemical reactions in the structure of DNA and RNA. Adjacent pyrimidine molecules form dimers and block the reproduction of bacteria, as a result, causing their death.

Destruction of microorganisms using UV radiation

The UV intensity needed for the elimination of microorganisms, such as yeasts, bacteria and viruses was previously investigated and reported by UVP Inc. A table below shows an amount of germicidal, shortwave (254 nm) UV energy needed for complete destruction of certain microorganisms.

Photochemical reaction

The diagram below shows the process of formation of pyrimidine dimers using thymine as an example (source: <http://www.photobiology.info>).



Bacteria organisms	Energy (mW/cm ² /s)
Bacillus anthracis	8.7
S. enteritidis	7.6
B. Megatherium sp. (veg.)	2.5
B. Megatherium sp. (spores)	5.2
B. paratyphosus	6.1
B. subtilis	11.0
B. subtilis spores	22.0

List continues on the next page ...

Other microorganisms	Energy (mW/cm ² /s)
Yeast	
Saccharomyces Ellipsoideus	13.2
Saccharomyces Sp.	17.6
Saccharomyces Cerevisiae	13.2
Brewer's Yeast	6.6
Baker's Yeast	8.8
Common Yeast Cake	13.2

List continues on the next page ...

Bacteria organisms

Energy (mW/cm²/s)

Clostridium tetani	22.0
Corynebacterium diphtheriae	6.5
Eberthella typosa	4.1
Escherichia coli	6.6
Micrococcus cadidus	12.3
Micrococcus sphaeroides	15.4
Mycobacterium tuberculosis	1.0
Neisseria catarrhalis	8.5
Phytomonas tumefaciens	8.5
Proteus vulgaris	6.6
Pseudomonas aeruginosa	10.5
Pseudomonas fluorescens	6.6
S. typhimuseum	15.2
Salmonella	10.0
Sarcina lutea	26.4
Serratia marcescens	6.1
Dysentery bacilli	4.2
Shigella paradysenteriae	3.2
Spirillum rubrum	6.1
Staphylococcus albus	5.7
Staphylococcus aureus	6.6
Streptococcus hemolyticus	5.5
Streptococcus lactis	8.8
Streptococcus viridans	3.8

Other microorganisms

Energy (mW/cm²/s)


Mold spores	
Penicillium Roqueforti	26.4
Penicillium Expansum	22.0
Penicillium Digitatum	88.0
Aspergillus Glaucus	88.0
Aspergillus Flavus	99.0
Aspergillus Niger	330.0
Rhisopus Nigricans	220.0
Mucor Racemosus A	35.2
Mucor Racemosus B	35.2
Oospora Lactis	11.0
Virus	
Bacteriophage (E. Coli)	6.6
Tobacco Mosaic	44.0
Influenza	6.6
Protozoa	
Paramecium	200.0
Nematode Eggs	92.0
Chlorella Vulgaris (Algae)	22.0

Table 1, Destruction chart of bacteria and various organisms (source: UVP Inc.)

UV Intensity measurements of Philips 25W bactericidal UV lamp


UV intensity depends on the distance from the UV source.

The graph below shows that UV intensity drops dramatically as the distance increases.



UV intensity, mW/cm ²	Distance, cm
18.0	0
9.3	2
5.0	5
2.8	10
2.2	15
1.7	20

Sensitivity of microorganisms to UV radiation intensity in UV air recirculators UVR-M and UVR-Mi



Yeast

Saccharomyces cerevisiae
Brewer's yeast
C. albicans
C. tropicalis
C. stellatoidea

Vegetative Bacteria

Clostridium tetani
Mycobacterium tuberculosis
Salmonella
Dysentery bacilli
Staphylococcus aureus
Streptococcus hemolyticus

Viruses

Bacteriophage (E. coli)
Influenza
Adenoviridae family
Retroviridae family
Coronaviridae family

Indoor pollution level before and after recirculator operation



ACKNOWLEDGEMENTS

This work was performed in R&D laboratory of company Biosan (Latvia, www.biosan.lv) and was supported by RSU Institute of Microbiology and Virology (Latvia).

Reproduction in any form is prohibited without the written permission of the administration of company Biosan SIA. Quoting is allowed only with reference to this source.

Creation: **September 2011**

Revision: **March 2020**