To SIA BIOSAN  
President  
Vasilijs Bankovskis

## Report

From 13.02.2018 to 13.10.2018, approbation of the following equipment produced by BioSan was conducted in SIA Kodolmedicīnas Klīnika radiopharmacy laboratory within the framework of the scientific project *Development and introduction in clinical practice of innovative methods of diagnostics and therapy of malignant tumours using molecularly targeted radionuclides produced in Latvia*: (1) laboratory centrifuge LMC-3000 with R12/10 rotor; (2) programmable mini-shaker Multi Bio 3D with PDM holder; (3) S-Bt Smart BioTherm compact CO2 incubator. This equipment has been used for standard in vitro test for the purpose of proving ability of mammal cancer cell lines to bind molecularly targeted radionuclides: S-Bt Smart BioTherm compact CO2 incubator has been used for keeping and incubating mammal cancer cell cultures, mini-shaker Multi Bio 3D for provision of mixing components of more uniform reaction, and laboratory centrifuge LMC-3000 for sedimentation and washing of cell suspension at various stages of reaction. All additions of reagents were performed by using BioSan automatic pipettes and disposable tips.

### *Specific circumstances of using the equipment:*

1. **LMC-3000**: Centrifuging has been running 3 minutes with 300 g or 900 g 15 ml Sarstedt in vials with screw caps. Speed provided by the equipment is suitable for cell experiments. It works perfectly, gathers speed fast, and also stops sufficiently fast. R12/10 rotor type suits needs of cell biology research.
2. CO2 Incubator **S-Bt Smart Biotherm**. Cells for experiments have been incubated for four or more hours at 37°C and in 5% CO2 atmosphere. Experiments used Sarstedt plastics suitable for mammal cell cultures flasks and 6 or 12 well plates. Equipment is operating well, keeps stable temperature within +/- 0.2°C. Chamber humidity during incubation has not been measured; water bath has been used during incubation.
3. **Biosan pipette** volume scope 0/5-10 ul, 2-20 ul, 20-200 ul, 100-1000 ul. Matched all our needs. Convenient in use, suitable tips.
4. **Multi Bio 3D** after adding the active substance, cells have been mixed in 6 or 12 well plate mode 10 rpm orbital rotation and 5° tilt change (vibro) mode. Using this equipment, we have conducted repeated in vitro test series three times (number of replicates for samples n=3), using five cancer cell lines. Obtained results have shown perfect replicability (degree of obtained data among individual replicates of analyzed objects was 94%) and reproducibility (experiments were conducted by different people at different times using the same equipment, match of data obtained in different experiments achieved 89%). No operation faults of used equipment have been observed in course of conducting experiments.

Member of the Board  
Andrejs Grīnbergs

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**EANM 18 / EP-0701**

Personalized medicine: method of in vitro determination of receptor binding for Lu-177 radio-pharmaceuticals therapy efficacy estimation

**Congress:** EANM 18

**Poster No.:** EP-0701

**Type:** Scientific e-Poster

**Keywords:** Clinical -> Therapy study -> Other treatments -> Local radionuclide treatment, 806 Local Radionuclide Treatment, 8. Radionuclide Therapy & Dosimetry

**Authors:** Z. Kalnina<sup>1</sup>, **T. Kusins**<sup>2</sup>, K. Svirks<sup>2</sup>, E. Rubena<sup>1</sup>, R. Kovaldins<sup>2</sup>, G. Kizane<sup>3</sup>, A. Berzina<sup>4</sup>, A. Grinbergs<sup>2</sup>; <sup>1</sup>Latvian Biomedical Research and Study centre, Riga, LATVIA, <sup>2</sup>Kodolmedicīnas klinika, Riga, LATVIA, <sup>3</sup>Institute of Chemical Physics, University of Latvia, Riga, LATVIA, <sup>4</sup>Riga East University Hospital, Clinical Department of Nuclear Medicine, Riga, LATVIA

**Aim**

Lu 177 -PSMA and Lu 177 -DOTA-TATE are currently used for the treatment of metastatic prostate cancer (mPC) and metastatic neuroendocrine tumours (mNET) when other therapeutic approaches fail, while their Ga 68 analogues help to precisely trace metastatic tumour spread. However, there is an interindividual heterogeneity of the targeted receptor expression levels even of a specific cancer type - e.g., among patients with mPC, in 5-10% there is no profound PSMA

...

**Authors:** Z. Kalnina<sup>1</sup>, **T. Kusins**<sup>2</sup>, K. Svirks<sup>2</sup>, E. Rubena<sup>1</sup>, R. Kovaldins<sup>2</sup>, G. Kizane<sup>3</sup>, A. Berzina<sup>4</sup>, A. Grinbergs<sup>2</sup>; <sup>1</sup>Latvian Biomedical Research and Study centre, Riga, LATVIA, <sup>2</sup>Kodolmedicīnas klinika, Riga, LATVIA, <sup>3</sup>Institute of Chemical Physics, University of Latvia, Riga, LATVIA, <sup>4</sup>Riga East University Hospital, Clinical Department of Nuclear Medicine, Riga, LATVIA

**i Disclosure 2 – Receiving Support (for details click on „i”)**

The research has been done in the framework of the project "Development and introduction of innovative methods in clinical practice for the diagnosis and treatment of malignant tumors using molecularly targeted radionuclides produced in Latvia", which is supported by philanthropist SIA "Mikrotikls" and administrated by the University of Latvia Foundation.

The conference-associated travel expenses of Zane Kalnina have been covered by a donation of Ekoosta, Ltd, administered by the University of Latvia Foundation.

For the project implementation, authors have received support in kind from a biotechnology company Biosan, Ltd and Arbor Medical Cooperation, Ltd.

Valdes loceklis  
Andrejs Grinbergs