

Steps

1

Incubation



● PSU-2T
*Mini-shaker for
 immunology*



● PST-60HL
*Plate Shaker-
 Thermostat*



● PST-60HL-4
*Plate Shaker-
 Thermostat*

2

Washing



● 3D-IW8
Inteliwasher



● FTA-2I
*Aspirator with Trap
 Flask*



ASSIST
pipette series



ASSISTBOY
pipette controller

New



ASSIST
*pipette series with
 volume lock*

3

Detection and analysis



● HIPO MPP-
 96
*Microplate
 Photometer*