

*LF-1, laboratory chest of drawers*



*Table T-4*



*Table T-4L-P*



*Table T-4L, with PCR-cabinet UVT-S-AR, 2 chests of drawers LF-1 and equipment for nucleic acid extraction*



*Table with shelves T-4L-P, chests of drawers LF-1 and equipment for nucleic acid extraction*



**Specifications T-4:**

Weight	23 kg
Overall dimensions (W×D×H)	800×600×770 mm

**Specifications T-4L:**

Weight	36 kg
Overall dimensions (W×D×H)	1290×600×770 mm

**Specifications T-4L-P:**

Weight	54 kg
Overall dimensions (W×D×H)	1290×600×1670 mm

**Specifications LF-1:**

Weight	28 kg
Overall dimensions (W×D×H)	300×450×705 mm

**Catalogue number:**

T-4	BS-040101-BK
T-4L	BS-040107-BK
T-4L-P	BS-050101-CK
LF-1	BS-050101-BK