



Use only tubes, plates suitable for centrifugation, adhere to the vessel manufacturer's specifications for optimal centrifugation parameters, including load and speed.



### STANDARD ROUND BOTTOM CENTRIFUGE TUBE

**10-15 ml**  
Ø: 16 mm  
L: 105 mm

×12 → **R-12/10**  
4,200 rpm  
BS-010208-BK

×24 → **R-24/10**  
4,000 rpm  
BS-010212-JK

### VACUTAINERS

**2-5 ml**  
Ø: 13 mm  
L: 82 mm  
[13 × 75]

×12 → **R-12/10**  
4,200 rpm  
BS-010208-BK

×24 → **R-24/10**  
4,000 rpm  
BS-010212-JK

BN-13/75

**4-8 ml**  
Ø: 13 mm  
L: 107 mm  
[13 × 100]

×12 → **R-12/10**  
4,200 rpm  
BS-010208-BK

×24 → **R-24/10**  
4,000 rpm  
BS-010212-JK

BN-13/100

**8-9 ml**  
Ø: 16 mm  
L: 105 mm  
[16 × 100]

×12 → **R-12/10**  
4,200 rpm  
BS-010208-BK

×24 → **R-24/10**  
4,000 rpm  
BS-010212-JK

BN-16/100

### MICROTUBES

**1.5-2 ml**  
Ø: 11 mm L: 40 mm

×24 → **RMT-24**  
6,000 rpm  
BS-010221-BK

### RACK FOR ROTORS

×1 → **RR-U**  
Rack for rotors  
BS-010208-UK

### STANDARD CONICAL BOTTOM CENTRIFUGE TUBE

**15 ml**  
Ø: 17 mm L: 120 mm

×12 → **R-12/15**  
4,200 rpm  
BS-010208-EK

**25 ml**  
Ø: 29 mm  
L: 77 mm

×6 → **R-6**  
4,200 rpm  
BS-010208-DK

×6 → **R-6P**  
4,200 rpm  
BS-010208-XK

B-25-6

**50 ml**  
Ø: 29 mm  
L: 115 mm

×6 → **R-6**  
4,200 rpm  
BS-010208-DK

×6 → **R-6P**  
4,200 rpm  
BS-010208-XK

### GEL CARDS

**Gel card**  
L: 74 mm  
H: 53 mm

×24 → **R-24GC**  
1,500 rpm  
BS-010208-VK

### 96, 384-WELL, DEEP-WELL PLATES, MICROPLATES

**Plate**  
L: 128 mm  
W: 85.6 mm  
H<sub>max</sub>: 45mm

×2 → **R-2**  
2,000 rpm  
BS-010208-AK

### SEMI- AND UNSKIRTED 96-WELL MICROPLATES

**Plate**  
L: 128 mm  
W: 85.6 mm  
H<sub>max</sub>: 45mm

×2 → **R-2**  
2,000 rpm  
BS-010208-AK

AP-96  
BS-010219-DK

### UNSKIRTED 384-WELL MICROPLATES

**Plate**  
L: 128 mm  
W: 85.6 mm  
H<sub>max</sub>: 45mm

×2 → **R-2**  
2,000 rpm  
BS-010208-AK

AP-384  
BS-010219-EK