



Use only tubes, plates suitable for centrifugation, adhere to the vessel manufacturer's specifications for optimal centrifugation parameters, including load and speed.



STANDARD ROUND BOTTOM CENTRIFUGE TUBE

10–15 ml
Ø: 16 mm
L: 105 mm

- Ri-12/10**
4,200 rpm
BS-010221-MK
- Ri-24/10**
4,000 rpm
BS-010221-LK
- BR-4U + BA-4/RT**
4,200 rpm
BS-010221-AK ⁺
BS-010221-3K ⁺

VACUTAINERS

2–5 ml
Ø: 13 mm
L: 82 mm
[13 × 75]

- Ri-12/10**
4,200 rpm
BS-010221-MK
- Ri-24/10**
4,000 rpm
BS-010221-LK
- BR-4U + BA-4/VT**
4,200 rpm
BS-010221-AK ⁺
BS-010221-4K ⁺

4–8 ml
Ø: 13 mm
L: 107 mm
[13 × 100]

- Ri-12/10**
4,200 rpm
BS-010221-MK
- Ri-24/10**
4,000 rpm
BS-010221-LK
- BR-4U + BA-4/VT**
4,200 rpm
BS-010221-AK ⁺
BS-010221-4K ⁺

8–9 ml
Ø: 16 mm
L: 105 mm
[16 × 100]

- Ri-12/10**
4,200 rpm
BS-010221-MK
- Ri-24/10**
4,000 rpm
BS-010221-LK
- BR-4U + BA-4/VT**
4,200 rpm
BS-010221-AK ⁺
BS-010221-4K ⁺

STANDARD CONICAL BOTTOM CENTRIFUGE TUBE

15 ml
Ø: 17 mm
L: 120 mm

- Ri-12/15**
4,200 rpm
BS-010221-KK
- BR-4U + BA-4/FT**
4,200 rpm
BS-010221-AK ⁺
BS-010221-YK ⁺

25 ml
Ø: 29 mm
L: 77 mm

- Ri-6**
4,200 rpm
BS-010221-HK
- Ri-6P**
4,200 rpm
BS-010221-IK
- BR-4U + BA-2/50**
4,200 rpm
BS-010221-AK ⁺
BS-010221-CK ⁺

50 ml
Ø: 29 mm
L: 115 mm

- Ri-6**
4,200 rpm
BS-010221-HK
- Ri-6P**
4,200 rpm
BS-010221-IK
- BR-4U + BA-2/50**
4,200 rpm
BS-010221-AK ⁺
BS-010221-CK ⁺

MICROTUBES

1.5–2 ml
Ø: 11 mm
L: 40 mm

- RMT-24**
6,000 rpm
BS-010221-BK

5 ml
Ø: 16 mm
L: 60 mm

- BR-4U + BA-14/2U**
4,200 rpm
BS-010221-AK ⁺
BS-010221-JK ⁺

5 ml
Ø: 16 mm
L: 60 mm


- BR-4U + BA-8/5**
4,200 rpm
BS-010221-AK ⁺
BS-010221-FK ⁺



Use only tubes, plates suitable for centrifugation, adhere to the vessel manufacturer's specifications for optimal centrifugation parameters, including load and speed.



96, 384-WELL, DEEP-WELL PLATES, MICROPLATES



→ (x2) →




Ri-2
2,000 rpm
BS-010221-GK


Plate dimension

L: 128 mm W: 85.6 mm H_{max}: 45mm

SEMI- AND UNSKIRTED 96-WELL MICROPLATES



→ (x2) →




Ri-2
2,000 rpm
BS-010221-GK

Plate dimension

L: 128 mm W: 85.6 mm H_{max}: 45mm

UNSKIRTED 384-WELL MICROPLATES



→ (x2) →




Ri-2
2,000 rpm
BS-010221-GK

Plate dimension

L: 128 mm W: 85.6 mm H_{max}: 45mm

GEL CARDS



→ (x24) →



Ri-24GC
1,500 rpm
BS-010221-NK

Dimension

L: 74 mm H: 53 mm

RACK FOR ROTORS



→ (x1) →



RR-U
Rack for rotors
BS-010208-UK