

## T-4, Table

### DESCRIPTION

New modular design of laboratory furniture provides flexibility and ease of use.

### SPECIFICATIONS


Maximum load	50 kg
Overall dimensions (W×D×H)	800 × 600 × 770 mm
Weight	23 kg



### CAT. NUMBER

BS-040101-BK	T-4, Table
--------------	------------

### ACCESSORIES



#### UVC/T-AR

BS-040102-AAA

DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box  
**UVC/T-AR** is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, plexiglas walls and ...

[read more](#)




#### UVC/T-M-AR

With built in socket and inlet for power cords  
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box  
**UVC/T-M-AR** is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, glass walls, working ...

[read more](#)



#### UVT-B-AR

With built in socket and inlet for power cords  
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box  
**UVT-B-AR** is designed for clean operations with DNA samples. UV-Cabinet provide protection against contamination.

Model is a bench-top type, made of metal framework, working surface made of ...

[read more](#)



#### LF-1

BS-050101-BK  
Drawer unit

New modular design of laboratory furniture provides flexibility and ease of use.