

## P-02/05, platform

### DESCRIPTION

For 24 tubes 0.5 ml and 48 tubes 0.2 ml.

### SPECIFICATIONS

Capacity

24 and 48

Microtubes volume

0.5 ml and 0.2 ml



### CAT. NUMBER

BS-010216-DK

Platform P-02/05



#### MPS-1

With built-in universal platform  
High-Speed Multi Plate Shaker

High-Speed Multi Plate Shaker **MPS-1** can be used in virtually any application by providing adjustable mixing of reagents in microtest plates, PCR plates, deepwell plates and test tubes (shaking tubes ...

[read more](#)