

P-5, shelf for pipettes

DESCRIPTION

P-5, shelf for pipettes


SPECIFICATIONS

Overall dimensions (W×D×H)

230 × 140 × 50 mm

Capacity

5



CAT. NUMBER

BS-040104-CK

P-5

ACCESSORIES



UVC/T-M-AR

With built in socket and inlet for power cords

DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box **UVC/T-M-AR** is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, glass walls, working ...

[read more](#)



UVT-B-AR

With built in socket and inlet for power cords
 DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box **UVT-B-AR** is designed for clean operations with DNA samples. UV-Cabinet provide protection against contamination.

Model is a bench-top type, made of metal framework, working surface made of ...

[read more](#)



UVT-S-AR

BS-040107-AAA
 DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box **UVT-S-AR** is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, glass walls, working ...

[read more](#)